

Aesculap Surgical Instruments

NOIR Supreme Cut Scissors



The next step in precision performance

Noir

SUPREME CUT SCISSORS


Introducing the next step in precision performance.

With Aesculap NOIR Supreme Cut Scissors you will experience a new dimension in precise cutting and dissecting performance. The Microtome cutting edge technology permits extremely fine and smooth cutting with minimal force. The cutting blade is precision ground to the sharpest angle.

The combination of the Durotip® tungsten carbide technology with the Aesculap low-reflection NOIR coating, results in an extremely high surface hardness and therefore high-wear resistance.

The Aesculap surface smoothing technology provides a low-wear, internal blade surface which results in increased blade life. The combined Aesculap technologies offer precision performance every time you use your NOIR Supreme Cut Scissors.

- Data matrix code technologies - allows fast machine-readable instrument tracking with instrument management software such as instacount® PLUS.



The extremely hard (2000HV) low-reflection NOIR coating enhances the wear resistance of the instruments and ensures a long life span, while providing high resistance against corrosion.



Precision Performance Every Time

Precision Performance Means *Razor Sharp*

- Microtome technology
 - Precision ground for the sharpest cutting edge


Every Time Means *Consistent High Performance*

- Durotip® TC technology
 - Tungsten carbide blades
- Smoothing technology
 - Low-wear internal blade surface
- Low-reflection NOIR coating
 - Superior wear resistance over full length



Noir

SUPREME CUT SCISSORS



Microtome technology ensures ultimate cutting edge precision for both smooth and wavecut blade options. The razor-sharp microtome blade edge offers unparalleled sharpness with minimal use of force.

Delicate Delicate



BC210B
straight
4-3/8"

BC211B
curved
4-3/8"



NOIR Iris
Dissecting Scissors



BC252B
straight
6-3/4"

BC253B
curved
6-3/4"



NOIR MAYO
Dissecting Scissors

Delicate



BC283B
curved
6-1/2"

BC284B
curved
6-1/2"



NOIR MAYO-LEXER
Dissecting Scissors



BC257B
curved
4-1/2"

BC261B
curved
5-3/4"

BC263B
curved
7"

BC271B
curved
7"

BC265B
curved
8"

BC275B
curved
8"


BC277B
curved
9"

NOIR Strabismus
Dissecting Scissors

NOIR METZENBAUM
Dissecting Scissors


Noir

SUPREME CUT SCISSORS



The particularly delicate shape permits smooth dissection even in very delicate tissue structures.

	Delicate			Delicate	Delicate
					
	1/1	1/1	1/1	1/1	1/1
	BC273B curved 7"	BC259B curved 5-3/4"	BC267B curved 9"	BC279B curved 10-1/4"	BC281B curved 11"
					
NOIR TOENNIS-ADSON Dissecting Scissors		NOIR BABY-METZENBAUM Dissecting Scissors	NOIR NELSON-METZENBAUM Dissecting Scissors		




Smooth Blade Option

Item No.	Description
BC210B	NOIR Iris Dissecting Scissors, delicate, straight, 4-3/8"
BC211B	NOIR Iris Dissecting Scissors, delicate, curved, 4-3/8"
BC252B	NOIR MAYO Dissecting Scissors, straight, 6-3/4"
BC253B	NOIR MAYO Dissecting Scissors, curved, 6-3/4"
BC257B	NOIR Strabismus Dissecting Scissors, curved, 4-1/2"
BC259B	NOIR BABY-METZENBAUM Dissecting Scissors, curved, 5-3/4"
BC261B	NOIR METZENBAUM Dissecting Scissors, curved, 5-3/4"
BC263B	NOIR METZENBAUM Dissecting Scissors, curved, 7"
BC265B	NOIR METZENBAUM Dissecting Scissors, curved, 8"
BC267B	NOIR NELSON-METZENBAUM Dissecting Scissors, curved, 9"
BC271B	NOIR METZENBAUM Dissecting Scissors, delicate, curved, 7"
BC273B	NOIR TOENNIS-ADSON Dissecting Scissors, delicate, curved, 7"
BC275B	NOIR METZENBAUM Dissecting Scissors, delicate, curved, 8"
BC277B	NOIR METZENBAUM Dissecting Scissors, delicate, curved, 9"
BC279B	NOIR NELSON-METZENBAUM Dissecting Scissors, delicate, curved, 10-1/4"
BC281B	NOIR NELSON-METZENBAUM Dissecting Scissors, delicate, curved, 11"
BC283B	NOIR MAYO-LEXER Dissecting Scissors, delicate, curved, 6-1/2"
BC284B	NOIR MAYO-LEXER Dissecting Scissors, curved, 6-1/2"

Noir

SUPREME CUT SCISSORS



W

Aesculap Wavecut technology prevents tissue from slipping out during the cutting process.

Delicate



BC211WB
curved
4-3/8"



NOIR Iris
Dissecting Scissors,
with Wavecut



BC253WB
curved
6-3/4"



NOIR MAYO
Dissecting Scissors,
with Wavecut

Delicate




BC283WB
curved
6-1/2"



NOIR MAYO-LEXER
Dissecting Scissors,
with Wavecut



BC284WB
curved
6-1/2"



BC257WB
curved
4-1/2"

NOIR Strabismus
Dissecting Scissors,
with Wavecut




BC261WB **BC263WB** **BC271WB** **BC265WB** **BC275WB** **BC277WB**
curved curved curved curved curved curved
5-3/4" 7" 7" 8" 8" 9"

NOIR METZENBAUM
Dissecting Scissors,
with Wavecut


Noir

SUPREME CUT SCISSORS



The particularly delicate shape permits smooth dissection even in very delicate tissue structures.

	Delicate			Delicate	Delicate
	1/1	1/1		1/1	1/1
BC273WB curved 7"		BC259WB curved 5-3/4"		BC267WB curved 9"	BC279WB curved 10-1/4"
				BC281WB curved 11"	
NOIR TOENNIS-ADSON Dissecting Scissors, with Wavecut		NOIR BABY-METZENBAUM Dissecting Scissors, with Wavecut		NOIR NELSON-METZENBAUM Dissecting Scissors, with Wavecut	



Wavecut Blade Option

Item No.	Description
BC211WB	NOIR Iris Dissecting Scissors, delicate, with Wavecut, curved, 4-3/8"
BC253WB	NOIR MAYO Dissecting Scissors, with Wavecut, curved, 6-3/4"
BC257WB	NOIR Strabismus Dissecting Scissors, with Wavecut, curved, 4-1/2"
BC259WB	NOIR BABY-METZENBAUM Dissecting Scissors, with Wavecut, curved, 5-3/4"
BC261WB	NOIR METZENBAUM Dissecting Scissors, with Wavecut, curved, 5-3/4"
BC263WB	NOIR METZENBAUM Dissecting Scissors, with Wavecut, curved, 7"
BC265WB	NOIR METZENBAUM Dissecting Scissors, with Wavecut, curved, 8"
BC267WB	NOIR NELSON-METZENBAUM Dissecting Scissors, with Wavecut, curved, 9"
BC271WB	NOIR METZENBAUM Dissecting Scissors, delicate, with Wavecut, curved, 7"
BC273WB	NOIR TOENNIS-ADSON Dissecting Scissors, delicate, with Wavecut, curved, 7"
BC275WB	NOIR METZENBAUM Dissecting Scissors, delicate, with Wavecut, curved, 8"
BC277WB	NOIR METZENBAUM Dissecting Scissors, delicate, with Wavecut, curved, 9"
BC279WB	NOIR NELSON-METZENBAUM Dissecting Scissors, delicate, with Wavecut, curved, 10-1/4"
BC281WB	NOIR NELSON-METZENBAUM Dissecting Scissors, delicate, with Wavecut, curved, 11"
BC283WB	NOIR MAYO-LEXER Dissecting Scissors, delicate, with Wavecut, curved, 6-1/2"
BC284WB	NOIR MAYO-LEXER Dissecting Scissors, with Wavecut, curved, 6-1/2"



AESCULAP®

All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogs and price lists and to take legal actions.

©2008 AESCULAP. ALL RIGHTS RESERVED. PRINTED IN THE USA. DOC637 Rev. A 4M 7/08
Aesculap is an equal opportunity employer

B | BRAUN
SHARING EXPERTISE

Aesculap, Inc.

3773 Corporate Parkway
Center Valley, PA 18034
Phone: 800-282-9000

www.aesculapusa.com